

TRINON
TITANIUM



Q-BONE-GRAFTING-SET

**Suitable for pre-implantological augmentation methods.
Such as onlay-plastic and
meshsupported bone-reconstruction**

a complete set includes:

Art.-No.	Description
10567	complete set including Steribox, all instruments and 70 pcs. of screws (without mesh)

All positions are available separate :

SCD1307	fixingscrew \varnothing 1,3 x 7 mm (1 Pck. = 10 pcs.)
SCD1309	fixingscrew \varnothing 1,3 x 9 mm (1 Pck. = 10 pcs.)
SCD1311	fixingscrew \varnothing 1,3 x 11 mm (1 Pck. = 10 pcs.)
SCD1313	fixingscrew \varnothing 1,3 x 13 mm (1 Pck. = 10 pcs.)
SCD1003	fixingscrew \varnothing 1,0 x 3 mm (1 Pck. = 10 pcs.)
SCD1005	fixingscrew \varnothing 1,0 x 5 mm (1 Pck. = 10 pcs.)
SCD1007	fixingscrew \varnothing 1,0 x 7 mm (1 Pck. = 10 pcs.)
10569	drill for fixingscrew \varnothing 1,3 mm 1 pcs.
10563	drill for fixingscrew \varnothing 1,0 mm 1 pcs.
10564	drill for boneblock 1 pcs.
10558	screwdriver endpiece \varnothing 1,3 mm 1 pcs.
10565	screwdriver endpiece \varnothing 1,0 mm 1 pcs.
10559	screwdriver handpiece 1 pcs.